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
Place the tube into one of the tube holder holes (on the side of the box).

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
Without touching the edges of the tube, add the buffer solution.

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
Insert the nasal swab into the right nostril (about 2 cm) and twist it 5 times (about 15 seconds). With the same swab, repeat the procedure in the left nostril.

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
Immediately after collection, insert the swab into the tube by rotating it in the solution 5 to 10 times while squeezing the tube.

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
Remove the swab by wringing it out to extract as much liquid as possible. Place the swab in its original packaging.

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
Only one line appears in the control area (C). No line appears in the test area (T). The sample does not contain SARS-CoV-2 virus. If the symptoms persist, repeat a test within a few days. If symptoms worsen, or you have any medical questions, contact your doctor. Respect social distancing guidelines.

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
If the control line (C) does not appear, the test result is not valid. Insufficient sample volume or an incorrect procedure are the most common reasons for missing control lines. Repeat the procedure using a new test or contact your doctor for any other diagnosis or medical question about you.

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Two lines appear, one in the control area (C) and another in the test area (T). The sample contains the SARS-CoV-2 virus. Isolate yourself and confirm the result as soon as possible with a PCR test.

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Associate your result with one of the following options :

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
Wait a full 15 minutes before reading the result. Do not wait more than 20 minutes.

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
Remove the test from its foil pouch and place it flat on a clean surface.

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Put 3 drops of the mixture in the S/R well.

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EN English

TEST SUMMARY

Ninonasal® self-test is a rapid, single-use immunochromatographic test for use as a diagnostic aid in the detection of SARS-CoV-2 antigen in a nasal sample. Ninonasal® self-test is intened for use by lay users in a private setting as a self-diagnosis.

TEST PRINCIPLE

Ninonasal® self-test includes a paper test strip inside a device. The test is performed by collecting a nasal sample extracted in a buffer solution, then placing 3 drops of this mixture on the test strip. When the test is performed correctly, up to 2 lines may appear on the paper. The control line (C) is used to confirm that the test has been used correctly and be visible if the migration has taken place. The test line (T) will only be visible if the applied sample contains the SARS-CoV-2 virus (responsible for COVID-19).

PERFORMING AND INTERPRETING

See detailed procedure on the back.

USE LIMITATIONS

- The test will produce negative results in the following conditions: the amount of SARS-CoV-2 virus in the sample is below the minimum detection limit of the test or the nasal sample was not collected correctly.

- Do not use nasal spray within 12 hours before performing the test.
- To have accurate results, do not use viscous samples or samples that contains visible blood.
- A negative result never rules out the possibility of infection with COVID-19. If the test result is negative and clinical symptoms persist, additional testing or other type of diagnosis should be performed.
- Failure to follow the instructions for use can lead to false positive or negative results.
- Infectivity cannot be ruled out on the sole basis of a negative test result.
- Positive or negative results must be confirmed by a healthcare professional.
- This test cannot determine the quantitative value of SARS-CoV-2 virus in the sample.
- The concentration of SARS-CoV-2 virus in the sample may vary during the course of the illness and may decrease below the detection limit of the test (especially when the duration of symptoms increases). Samples collected within 7 days of onset of the first symptoms are generally more likely to give a positive result.

- The test confirms the presence of SARS-CoV-2 core antigen (N protein) in a symptomatic patient. If the result is negative, it should not be used as the sole criterion for diagnosing a SARS-CoV-2 infection, consult your doctor.
- The user should not make any medically important decision regarding his health state without first consulting a healthcare professional.
- The test is intended for adults. If the test is performed on children under 18 years old, it must be performed under the supervision of an adult.

WARNINGS AND PRECAUTIONS

- The test is for single use only. The test device, swab, tube as cap should not be reused.
- Do not use if the expiry date has passed.
- Do not use if the foil pouch is damaged (eg. torn, punctured) or opened.
- Do not open the aluminum pouch before reading the instructions carefully and being ready to test.
- Do not use any solution other than the buffer solution that was supplied with the test. Do not use an open or leaking buffer solution in bottle or unit-dose.
- Avoid eye/skin contact with the buffer solution.
- The test should be placed on a flat and clean surface while awaiting the result. The test should never face upwards.
- Do not eat, drink or smoke in the test handling area.
- If the instructions are not followed correctly, the results may be false.
- In case the pouch has been stored at 4-8°C, wait

at least 10 minutes for the test to return to room temperature (15-25°C).

STORAGE

- The test should be stored between 4°C and 30°C.
- Do not freeze.
- Do not open the foil pouch until you are ready to perform the test. The test should be performed immediately after opening.

DISPOSAL

After use, place the test and the kit content in its original package and dispose it according to local regulations.

PERFORMANCES

The sensitivity (or ability to detect true positives) of the test is 98% (confidence interval: IC95% [93-100]). The specificity (or ability to detect true negatives) of the test is 99% (confidence interval: IC95% [98-100]). The test results were compared with the results obtained with the standard method (RT-PCR).

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